

ultra-trace mercury analysis



Tekran Inc.

NEW DIRECTIONS IN ENVIRONMENTAL TECHNOLOGY

Series 2600



The most sensitive mercury system available today.



The Tekran® Series 2600 offers unparalleled performance and flexibility for those who require truly sensitive total mercury analysis.

The system can implement a wide range of different analytical techniques, depending on the requirements of your particular application. Unlike conventional systems, the **Series 2600** is available in a wide range of configurations, with options to suit your budget. The system can perform virtually any type of ultra-trace total mercury analysis. The **Series 2600** allows easy migration from a starter system to a fully automated, high throughput configuration.

Highlights

- Sensitivity measured in *parts per quadrillion*. Typical MDL: $<0.05 \text{ ppt}^1$
- Atomic Fluorescence based: Greater sensitivity, selectivity and dynamic range than AA based systems
- Unprecedented modularity and flexibility

Typical Applications

Total Mercury in Liquids

- Automated analysis of aqueous samples
 - *Dual stage gold preconcentration*
 - *Single stage gold preconcentration*
 - *Direct reading*
- May be used with stannous chloride or sodium borohydride reductant

Mercury in Gases


- Automated analysis of adsorbent sample cartridges using dual stage pre-concentration with thermal desorption

Typical Full Flow Diagram

Automated liquid analysis using dual stage gold pre-concentration.


(US EPA Method 1631)

¹ Dual stage pre-concentration. Analysis conducted in clean room with controlled blanks. Sample aliquot: 30 ml.



METHODS

Liquid Analysis via Dual Stage Gold Preconcentration US EPA Method 1631




The combination of dual stage preconcentration and atomic fluorescence detection provides the most sensitive, selective and resilient mercury analysis available today. After passing through the gas/liquid phase separator, the carrier gas is passed through a desiccant dryer and on to a gold coated sample collection cartridge.

This cartridge is then thermally desorbed onto an analytical gold trap. This


dual stage preconcentration method eliminates virtually all interfering compounds, ensuring excellent analytical results.

The solenoid valve **V2** ensures precise loading of the first cartridge. Valve **V1** vents effluent during loading, eliminating contamination of downstream components.

Peak obtained from 0.5 ppt standard using gold preconcentration



Liquid Analysis via Direct Atomic Fluorescence US EPA Method 245.7




The direct measurement technique provides good performance for many applications and allows for a faster analytical cycle. Carrier gas from the gas/liquid phase separator is passed through a Perma Pure® diffusion dryer and sent directly to the AF detector.

Valve **V2** provides precise introduction of the loaded carrier gas into the


detector and eliminates the need to use solenoids for switching liquid sample.

Systems equipped for dual stage preconcentration can be converted to and from the direct method in only minutes!

Five ppt standard using direct method



Dual Stage Gas Phase Analysis




The **Model 2600**, without any of the liquid handling components, can be used for gas analysis. Sample cartridges are inserted and the entire thermal desorption process is performed automatically.

The arrangement shown is optimized for gas analysis and allows the field cartridge to be changed while the analytical trap is still being analysed.

The solenoids vent cartridge effluent during preheating to eliminate downstream contamination.

Peak from gas sample with 1.7 ng/m³ of mercury (50 L integrated sample)



MODULARITY

MODEL 2600

Preconcentrator / Detector

The Tekran® Model 2600 is the most sensitive and flexible atomic fluorescence mercury analytical unit available today. The wide range of options available allows a system to be configured to suit virtually any application and budget.

ULTRA-SENSITIVE ATOMIC FLUORESCENCE DETECTOR

- Based on the Tekran tried and proven Cold Vapor Atomic Fluorescence Spectrophotometer (CVAFS) elemental mercury detector
- Very high sensitivity, (MDL < 0.1 pg)
- Wide dynamic range (>10⁵)
- Inherently linear
- Superior selectivity

FULL FEATURED DUAL STAGE PRECONCENTRATOR

- Individual cooling fans for rapid cycle times
- Continuously variable heater power prevents condensation and allows precleaning of cartridges
- Heating stations can be populated with virtually any type of preconcentration cartridge:
 - Pure gold, gold coated sand, coated silica-gel, coated glass beads, etc.

MULTIPLE SCRUBBER / DRYER OPTIONS

- Desiccant dryer/acid gas scrubber cartridge (eg: sodalime)
- Perma Pure® diffusion dryer

FLOW CONTROL SOLENOID VALVES

- Multi-function all Teflon solenoid valves:
 - Provide precise analytical loadings
 - Ensure clean sample is sent to the detector by venting flows during preclean operations
 - Separate flow paths to allow overlapped adsorption and desorption operations
- Only the gas phase is switched by the solenoids
 - Eliminates contamination and carry over problems inherent in switching liquid sample streams
 - Prevents valve damage caused by crystallization of reductant within valve body

HI-RESOLUTION DATA CAPTURE

- 24 bit A/D convertor ensures accurate measurement of both large and small peaks

FLEXIBLE FLOW PATH CONFIGURATION

- All flow path components are made of Teflon®.
- Convenient push-on fittings allow easy component replacement and flow path changes.

MASS FLOW CONTROLLER

- Optional precision mass flow controller may be used to set critical analytical flows
- Flow rates may be automatically programmed over the course of each analysis

BALL FLOW CONTROLLERS

- Up to three ball flowmeters with manual adjustment valves may be configured on the front panel
- Interchangeable tube sets allow full scale flows of 50, 100, 200, or 500 ml/min, 1, 2 or 5 l/min

COMMUNICATIONS

- All system functions (including the autosampler) are performed through a single serial port
- A notebook PC may serve as the controller for portable applications

EXTERNAL EVENT CONTROL

- Two opto-isolated external outputs can control external equipment

HIGH EFFICIENCY PHASE SEPARATOR

- Removable center rod for easy cleaning
- Low dead volume
- High efficiency



Low cost version for use with external Model 2500 detector



Version with internal detector and two ball flowmeters. (no MFC)



Full configuration with internal AF detector, mass flow controller, and two ball flowmeters





MODEL 2620

Automatic Sample Changer

The Model 2620 Automatic Sample Changer (Auto Sampler), is a practical, cost effective solution for laboratories with large numbers of samples to run. The addition of our Tekran-MDS data handling system creates a fully automated, high throughput mercury analysis system. For the first time, users have the ability to accurately perform ultra-trace mercury analysis in the parts per quadrillion range, without user intervention.

EASY INTEGRATION

- Connects directly to the Model 2600
- Does not require a second serial port on your PC

WASH STATION

- Glass recirculating wash station provides continuous rinsing of sample probe
- Provides clean deionized water for system flushing between samples
- May be operated in either recirculating or single pass mode

RACKS

- Auto sampler holds up to four individual rack modules
- Usable with a wide range of industry standard rack types
- Different rack types may be mixed within a single run
- Two typical rack types are shown below.

Rack code	Positions per rack	Total positions	Volume (ml)	Material
22	44	176	30	Borosilicate glass
112	12	48	60	Teflon



Control Software

The Tekran Mercury Data System (Tekran-MDS) is a full featured Windows® based application that provides all control, display and data reduction functions for the system.

SIMPLE DATA ENTRY

- Samples and standards are entered into a run worksheet. Detailed information for each sample is available in a single sample spreadsheet.

FULLY PROGRAMMABLE SAMPLE CYCLE


- Convenient editor allows full user control of analytical cycle, including:
 - Pump speed, carrier flow rate, wash pump, solenoid activations, heater power levels, cooling fans, peak acquisition, autosampler operations

PEAK INTEGRATION

- Detector output for each run is automatically captured and integrated
- Manual reintegration is possible using graphical user interface

CALIBRATION

- Multi-point least squares fit
- Variety of curve fitting options
- Export of analytical results to your own QA/QC packages or LIMS



REPORTING

- Report formats may be user configured
- Easy export to other applications

MODEL 2610

Peristaltic Pump Module

The addition of the Model 2610 Pump Module creates a complete atomic fluorescence analyzer for mercury in aqueous samples. This combination allows you to analyze background levels in virtually any type of liquid sample. The Pump Module fits directly beneath the Model 2600 and the two together require only 9" of bench space.

WASH PUMP


- Provides continuous source of DI water for the wash station on the Automatic Sample Changer
- Two channel peristaltic, fixed speed
- Automatic shutoff after run is complete

MAIN PUMP

- Five channel variable speed precision peristaltic pump
- Precise flow delivery and long term stability is ensured by using closed loop optical speed controller
- Pump speed may be set either manually or automatically under program control
- Variable speed control minimizes use of reagents and sample while allowing quick analysis times

Series 2600

SPECIFICATIONS

MODEL 2600	Control Unit	MODEL 2610	Pump Module	MODEL 2620	Automatic Sampler
DETECTOR Type: Atomic fluorescence, 253.7 nm Sensitivity: < 0.1 pg Warm up Time: < 10 min PRECONCENTRATOR Stages: 2 max. FLOW REGULATION Mass Flow Meter: 1 max. Ball Flow Meter: 3 max. DIMENSIONS Width: 9" (230 mm) Height: 12" (300 mm) Depth: 17" (430 mm) WEIGHT Instrument: 30 lb (14 kg) Shipping: 45 lb (20 kg) POWER REQUIREMENTS 100-120, 220-240 VAC, 50-60 Hz., 200 VA max.		MAIN PUMP Type: 5 channel peristaltic Speed: 0 - 200 rpm Regulation: closed loop optical tachometer Control: local or computer (continuously variable) WASH PUMP Type: 2 channel peristaltic Speed: 100 rpm, fixed Control: local or computer DIMENSIONS Width: 9" (230 mm) Height: 5" (130 mm) Depth: 17" (430 mm) WEIGHT Instrument: 14 lb (7 kg) Shipping: 24 lb (11 kg) POWER REQUIREMENTS 100-120, 220-240 VAC, 50-60 Hz., 100 VA max.		SAMPLER UNIT Type: XYZ 3 axis with four stationary rack modules Capacity: Code 112 - 48 positions (60 ml) Code 22 - 176 positions (30 ml) (Other rack and tube types available) Control: Serial RS-232. Probe Rinse: Dedicated recirculating DI water rinse station. Arm Speed: X: 25.15 cm/sec; Y: 25.25 cm/sec; Z: 2.0, 3.0, 6.2, 12.7, or 24.7 cm/sec POSITIONING Accuracy: +/- 1 mm in XYZ dimensions Repeatability: +/- 0.25 mm in XYZ dimensions WORKING AREA X: 13" (33 cm); Y: 9.5" (24.1 cm) DIMENSIONS Width: 21" (535 mm) Depth: 17" (430 mm) Height: 14.6" (371 mm) Height with arm: 22" (560 mm) WEIGHT Instrument: 38 lb (17 kg) Shipping: 49 lb (22 kg) POWER REQUIREMENTS 100-120, 220-240 VAC, 50-60 Hz., 75 VA max.	
Tekran® trademark of Tekran Inc. Windows® trademark of Microsoft Corp. Perma Pure® trademark of Perma Pure Inc. Due to continuing product development, specifications are subject to change without notice.					

For more information contact:

Tekran Inc.

132 Railside Road, Unit 1, Toronto, Ontario Canada M3A 1A3

Tel: +1 416 449-3084 Fax: +1 416 449-9298

Toll Free Sales Line (US and Canada):

1-888-3TEKRAN (1-888-383-5726)

AUTHORIZED REPRESENTATIVE:

For the latest products and developments, visit our website at: <http://tekran.com>

Tekran Inc.

info@tekran.com