

# Cole-Parmer Instrument Company

**ThermoFisher**

Manufactured by: **SCIENTIFIC**

## INFORMATION SHEET

### Flow Rate (mL/min) and Ordering Information 115V AC 60 Hz

Motor	Motor rpm	Current Rating (Ampere)	Pump Head						
			D-2533-0001	-0002	-0004	-0003	-0005	-0007	-0006
			13	14	16	15	17	24	18
B-3101-0001	1	0.25	0.1 mL/min 900-0809	0.2 mL/min 900-0816	0.8 mL/min 900-0823	1.6 mL/min 900-0830	2.8 mL/min 900-0837	2.8 mL/min 900-0844	3.8 mL/min 900-0851
B-3101-0002	6	0.15	0.4 900-0810	1.3 900-0817	4.8 900-0824	9.6 900-0831	16.8 900-0838	16.8 900-0845	22.8 900-0852
B-3101-0003	30	0.6	1.8 900-0811	5.3 900-0818	24 900-0825	48 900-0832	84 900-0839	84 900-0846	114 900-0853
B-3101-0004	80	0.6	4.8 900-0812	16.8 900-0819	64 900-0826	128 900-0833	224 900-0840	224 900-0847	304 900-0854
B-3101-0005	150	0.6	9 900-0813	31.5 900-0820	120 900-0827	240 900-0834	420 900-0841	420 900-0848	570 900-0855
B-3101-0006	200	0.6	12 900-0814	42 900-0821	160 900-0828	320 900-0835	560 900-0842	560 900-0849	760 900-0856
B-3101-0007	250	0.6	15 900-0815	52.5 900-0822	200 900-0829	400 900-0836	700 900-0843	700 900-0850	950 900-0857

→ PBM  
- cal

### Cole-Parmer

1-800-MASTERFLEX (627-8373) (U.S. and Canada only) • 11 (847) 549-7600 (Outside U.S.)  
(847) 549-7600 (Local) • [www.masterflex.com](http://www.masterflex.com) • [techinfo@coleparmer.com](mailto:techinfo@coleparmer.com)

### Thermo Fisher Scientific

1-800-637-3739 (U.S. and Canada only) • 11 (847) 381-7050 (Outside U.S.)  
847-381-7050 (Local) • [www.thermoscientific.com](http://www.thermoscientific.com) • [fluidhandling@thermoscientific.com](mailto:fluidhandling@thermoscientific.com)